

## **Metaldehyde (Homopolymer) 100 µg/mL in Acetonitrile**

Art. ID LM24-N-15900-5809-100AN5  
Unit 5 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1648 II E2

### Description

Conservation: 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Metaldehyde (Homopolyme r)	[9002-91-9]	100	µg/mL		