

Cyantraniliprole IN-JSE76 10 µg/mL in Acetonitrile

Art. ID LM24-N-15900-5851-10AN5
Unit 5 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1648 II E2

Description

Conservation: 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Cyantraniliprole IN-JSE 76	[n/a]	10	µg/mL		