

Folpet-d4 100 µg/mL in Acetonitrile

Art. ID LM24-N-15900-5857D-100AN1
Unit 1 mL
Deliverydetails De Minimis UN1648 II E2

Description

Conservation: 18 months (originally closed bottle).

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|----------------|---------------------|-------|--------|--------|
| | Folpet D4 | [1327204-12-5] | 100 | µg/mL | | |