

2,5-Dichlorobenzoic acid-methyl ester 100 µg/mL in Acetonitrile

Art. ID LM24-N-15900-5881-100AN5
Unit 5 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1648 II E2

Description

Conservation: 18 months (originally closed bottle).

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|---|-------------|---------------------|-------|--------|--------|
| | 2,5-Dichlorobenzoic aci d-methyl ester | [2905-69-3] | 100 | µg/mL | | |