

**(±)-MDMA-d5 10 µg/mL in Methanol**

Art. ID LM24-N-18000-7811D-10ME5  
Unit 5 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

**Description**

Conservation: 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	(±)-MDMA-D5	[136765-43-0]	10	µg/mL		