

(±)-MDMA-d5 10 µg/mL in Methanol

Art. ID LM24-N-18000-7811D-10ME5
Unit 5 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

Description

Conservation: 18 months (originally closed bottle).

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|---------------|---------------------|-------|--------|--------|
| | (±)-MDMA-D5 | [136765-43-0] | 10 | µg/mL | | |