

### **4-Formylaminoantipyrine-D3 10 µg/mL in Acetonitrile**

Art. ID LM24-N-18000-7875D-10AN5  
Unit 5 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1648 II E2

#### Description

Conservation: 18 months (originally closed bottle).

| Text/Information | Analyte/Parameter              | CAS number | Concentration/Value | Unit  | Method | Source |
|------------------|--------------------------------|------------|---------------------|-------|--------|--------|
|                  | 4-Formylaminoantipyrine<br>-d3 |            | 10                  | µg/mL |        |        |