

2,2'-Diaminodiphenylmethane 50 µg/mL in Acetonitrile

Art. ID LM24-N-19000-0465-50AN5
Unit 5 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1648 II E2

Description

Conservation: 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2,2'-Methylenedianiline	[6582-52-1]	50	µg/mL		