

2-Ethylhexanoic acid-d15 1000 µg/mL in Acetone

Art. ID LM24-N-20000-0149D-1000AC5
Unit 5 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1090 II E2

Description

Conservation: 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	(+/-)-2-Ethylhexanoic-D 15 Acid	[352431-38-0]	1000	µg/mL		