

2-Ethylhexanoic acid-d15 100 µg/mL in Acetone

Art. ID LM24-N-20000-0149D-100AC5
Unit 5 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1090 II E2

Description

Conservation: 18 months (originally closed bottle).

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|------------------------------------|---------------|---------------------|-------|--------|--------|
| | (+/-)-2-Ethylhexanoic-D 15 Acid | [352431-38-0] | 100 | µg/mL | | |