

## **2-Ethylhexanoic acid-d15 100 µg/mL in Acetone**

Art. ID LM24-N-20000-0149D-100AC5  
Unit 5 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1090 II E2

### Description

Conservation: 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	(+/-)-2-Ethylhexanoic-D 15 Acid	[352431-38-0]	100	µg/mL		