

Custom Solution - 9 components 100 mg/L each of in Methanol

Art. ID LM24-CUS-72764
Unit 1 x 1 mL ampoule
Deliverydetails De Minimis UN1230 (6.1) II E2

Description

Conservation 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	N-Nitroso-dimethylamine	[62-75-9]	100	mg/L		
	N-Nitroso-diethylamine	[55-18-5]	100	mg/L		
	N-Nitroso-di-n-propylamine	[621-64-7]	100	mg/L		
	N-Nitroso-di-n-butylamine	[924-16-3]	100	mg/L		
	N-Nitrosopiperidine	[100-75-4]	100	mg/L		
	N-Nitrosopyrrolidine	[930-55-2]	100	mg/L		
	N-Nitrosomorpholine	[59-89-2]	100	mg/L		
	N-Nitroso-N-methylaniline	[614-00-6]	100	mg/L		
	N-Nitroso-N-ethylaniline	[612-64-6]	100	mg/L		