

Custom Solution - 27 components in HNO3 5%

| | |
|-----------------|--|
| Art. ID | LM24-CUS-72765 |
| Unit | 1 x 500 mL bottle |
| Deliverydetails | ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2 |

Description

Conservation 24 months (originally closed bottle).

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Silver (Ag) | [7440-22-4] | 100 | mg/L | | |
| | Aluminium (Al) | [7429-90-5] | 100 | mg/L | | |
| | Arsenic (As) | [7440-38-2] | 100 | mg/L | | |
| | Boron (B) | [7440-42-8] | 1000 | mg/L | | |
| | Barium (Ba) | [7440-39-3] | 100 | mg/L | | |
| | Bismuth (Bi) | [7440-69-9] | 100 | mg/L | | |
| | Calcium (Ca) | [7440-70-2] | 2000 | mg/L | | |
| | Cadmium (Cd) | [7440-43-9] | 100 | mg/L | | |
| | Cobalt (Co) | [7440-48-4] | 100 | mg/L | | |
| | Chromium (Cr) | [7440-47-3] | 100 | mg/L | | |
| | Copper (Cu) | [7440-50-8] | 100 | mg/L | | |
| | Iron (Fe) | [7439-89-6] | 200 | mg/L | | |
| | Potassium (K) | [7440-09-7] | 1500 | mg/L | | |
| | Lithium (Li) | [7439-93-2] | 100 | mg/L | | |
| | Magnesium (Mg) | [7439-95-4] | 1000 | mg/L | | |
| | Manganese (Mn) | [7439-96-5] | 100 | mg/L | | |
| | Molybdenum (Mo) | [7439-98-7] | 100 | mg/L | | |
| | Sodium (Na) | [7440-23-5] | 1500 | mg/L | | |
| | Nickel (Ni) | [7440-02-0] | 100 | mg/L | | |
| | Lead (Pb) | [7439-92-1] | 100 | mg/L | | |
| | Antimony (Sb) | [7440-36-0] | 100 | mg/L | | |
| | Selenium (Se) | [7782-49-2] | 100 | mg/L | | |
| | Strontium (Sr) | [7440-24-6] | 100 | mg/L | | |
| | Titanium (Ti) | [7440-32-6] | 100 | mg/L | | |
| | Thallium (Tl) | [7440-28-0] | 100 | mg/L | | |
| | Vanadium (V) | [7440-62-2] | 100 | mg/L | | |
| | Zinc (Zn) | [7440-66-6] | 100 | mg/L | | |