

Custom Solution - 3 components in Water

Art. ID LM24-CUS-72767
Unit 1 x 500 mL bottle
Deliverydetails No Dangerous Good /not restricted

Description

Conservation 24 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Phosphorus (P)	[7723-14-0]	600	mg/L		
	Sulfur (S)	[7704-34-9]	1000	mg/L		
	Silicon (Si)	[7440-21-3]	1000	mg/L		