

Custom Solution - 12 components 100 mg/L each of in Methanol

Art. ID LM24-CUS-72768
Unit 3 x 1 mL ampoule
Deliverydetails De Minimis UN1230 (6.1) II E2

Description

Conservation 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Ethylacetate	[141-78-6]	100	mg/L		
	Benzylalcohol	[100-51-6]	100	mg/L		
	2-Butanone	[78-93-3]	100	mg/L		
	Cyclohexane	[110-82-7]	100	mg/L		
	Ethanol	[64-17-5]	100	mg/L		
	2-Propanol	[67-63-0]	100	mg/L		
	Acetic acid-methyl ester	[79-20-9]	100	mg/L		
	r					
	2-Methyl-1-propanol	[78-83-1]	100	mg/L		
	4-Methyl-2-pentanone	[108-10-1]	100	mg/L		
	n-Hexane	[110-54-3]	100	mg/L		
	n-Heptane	[142-82-5]	100	mg/L		
	Butylacetate	[123-86-4]	100	mg/L		