

Pharma - Mix 18 - 0.0000025 M each of in Water

Art. ID LM24-CUS-10293-2X10ML
Unit 2 x 10 mL
Deliverydetails No Dangerous Good /not restricted

Description

Conservation: 18 months (originally closed bottle). NOT ISO 17034 CERTIFIED.

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|------------|---------------------|------|--------|--------|
| | L-Alanine | [56-41-7] | 0,0000025 | M | | |
| | L-Arginine | [74-79-3] | 0,0000025 | M | | |
| | Aspartic acid | [56-84-8] | 0,0000025 | M | | |
| | L-Cystine | [56-89-3] | 0,0000025 | M | | |
| | L-Glutamic acid | [56-86-0] | 0,0000025 | M | | |
| | Glycine | [56-40-6] | 0,0000025 | M | | |
| | L-Histidine | [71-00-1] | 0,0000025 | M | | |
| | L-Isoleucine | [73-32-5] | 0,0000025 | M | | |
| | L-Leucine | [61-90-5] | 0,0000025 | M | | |
| | L-Lysine | [56-87-1] | 0,0000025 | M | | |
| | Methionine | [63-68-3] | 0,0000025 | M | | |
| | Phenylalanine | [63-91-2] | 0,0000025 | M | | |
| | L(-)-Proline | [147-85-3] | 0,0000025 | M | | |
| | L-Serine | [56-45-1] | 0,0000025 | M | | |
| | L-Threonine | [72-19-5] | 0,0000025 | M | | |
| | Tyrosine | [60-18-4] | 0,0000025 | M | | |
| | L-Tryptophan | [73-22-3] | 0,0000025 | M | | |
| | L-Valine | [72-18-4] | 0,0000025 | M | | |