

Austenitic Stainless Steel, UNS S32180, Chips

Art. ID MBH-13X32180A-C
Unit 50 g (chips)
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	0,043	%		
	Arsenic (As)	[7440-38-2]	~0,003	%		
	Boron (B)	[7440-42-8]	~0,0011	%		
	Carbon (C)	[7440-44-0]	0,031	%		
	Cobalt (Co)	[7440-48-4]	0,04	%		
	Chromium (Cr)	[7440-47-3]	18,92	%		
	Copper (Cu)	[7440-50-8]	0,49	%		
	Iron (Fe)	[7439-89-6]	Base	%		
	Manganese (Mn)	[7439-96-5]	2,11	%		
	Molybdenum (Mo)	[7439-98-7]	0,245	%		
	Nitrogen (N)	[7727-37-9]	0,0067	%		
	Niobium (Nb)	[7440-03-1]	~0,0021	%		
	Nickel (Ni)	[7440-02-0]	10,16	%		
	Phosphorus (P)	[7723-14-0]	0,007	%		
	Sulfur (S)	[7704-34-9]	0,0093	%		
	Antimony (Sb)	[7440-36-0]	~0,0011	%		
	Silicon (Si)	[7440-21-3]	0,485	%		
	Tin (Sn)	[7440-31-5]	0,0116	%		
	Titanium (Ti)	[7440-32-6]	0,81	%		
	Vanadium (V)	[7440-62-2]	0,026	%		
	Tungsten (W)	[7440-33-7]	0,039	%		