

Manganese Brass, Chips

Art. ID MBH-31XB25C-C
Unit 50 g (chips)
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Silver (Ag) | [7440-22-4] | 0,0025 | % | | |
| | Aluminium (Al) | [7429-90-5] | 0,48 | % | | |
| | Arsenic (As) | [7440-38-2] | 0,027 | % | | |
| | Bismuth (Bi) | [7440-69-9] | 0,061 | % | | |
| | Cadmium (Cd) | [7440-43-9] | 0,00019 | % | | |
| | Chromium (Cr) | [7440-47-3] | ~0,0002 | % | | |
| | Copper (Cu) | [7440-50-8] | 55,6 | % | | |
| | Iron (Fe) | [7439-89-6] | 0,014 | % | | |
| | Manganese (Mn) | [7439-96-5] | 0,128 | % | | |
| | Nickel (Ni) | [7440-02-0] | 0,255 | % | | |
| | Phosphorus (P) | [7723-14-0] | 0,116 | % | | |
| | Lead (Pb) | [7439-92-1] | 0,272 | % | | |
| | Antimony (Sb) | [7440-36-0] | 0,075 | % | | |
| | Silicon (Si) | [7440-21-3] | 0,235 | % | | |
| | Tin (Sn) | [7440-31-5] | 0,62 | % | | |
| | Zinc (Zn) | [7440-66-6] | 42,2 | % | | |