

Tin-base (Lead-Free) Solder, 50g Chips

Art. ID MBH-74XCA1-21-C
Unit 50 g (chips)
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	0,032	%		
	Antimony (Sb)	[7440-36-0]	0,019	%		
	Bismuth (Bi)	[7440-69-9]	0,014	%		
	Cadmium (Cd)	[7440-43-9]	0,006	%		
	Copper (Cu)	[7440-50-8]	0,681	%		
	Gold (Au)	[7440-57-5]	0,0047	%		
	Indium (In)	[7440-74-6]	0,0031	%		
	Iron (Fe)	[7439-89-6]	0,0026	%		
	Lead (Pb)	[7439-92-1]	0,079	%		
	Nickel (Ni)	[7440-02-0]	0,0003	%		
	Phosphorus (P)	[7723-14-0]	~34	ppm		
	Silver (Ag)	[7440-22-4]	0,44	%		
	Zinc (Zn)	[7440-66-6]	~6	ppm		