

Lead-Tin Solder Binary Alloy, Ø40x8mm

Art. ID MBH-91XS63P-22-G
Unit disc
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silver (Ag)	[7440-22-4]	0,0046 ± 0,0004	%		
	Bismuth (Bi)	[7440-69-9]	0,0048 ± 0,0008	%		
	Copper (Cu)	[7440-50-8]	0,0016 ± 0,0003	%		
	Iron (Fe)	[7439-89-6]	0,0020 ± 0,0004	%		
	Indium (In)	[7440-74-6]	0,0044 ± 0,0006	%		
	Nickel (Ni)	[7440-02-0]	0,0006 ± 0,0003	%		
	Antimony (Sb)	[7440-36-0]	0,018 ± 0,004	%		
	Tin (Sn)	[7440-31-5]	62,2 ± 0,2	%		
	Arsenic (As)	[7440-38-2]	~8	ppm		
	Gold (Au)	[7440-57-5]	~4	ppm		
	Calcium (Ca)	[7440-70-2]	~2410	ppm		
	Cadmium (Cd)	[7440-43-9]	~2	ppm		
	Mercury (Hg)	[7439-97-6]	~34	ppm		
	Magnesium (Mg)	[7439-95-4]	~6	ppm		
	Lead (Pb)	[7439-92-1]	Balance	ppm		
	Selenium (Se)	[7782-49-2]	~52	ppm		
	Tellurium (Te)	[13494-80-9]	~64	ppm		
	Zinc (Zn)	[7440-66-6]	~16	ppm		