

Tin-base (Lead-Free) Solder, 50g Chips

Art. ID MBH-74XCA5-21-C
Unit 50 g (chips)
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	0,0004	%		
	Antimony (Sb)	[7440-36-0]	0,133	%		
	Arsenic (As)	[7440-38-2]	0,011	%		
	Bismuth (Bi)	[7440-69-9]	0,0056	%		
	Cadmium (Cd)	[7440-43-9]	0,0019	%		
	Copper (Cu)	[7440-50-8]	1,22	%		
	Gold (Au)	[7440-57-5]	0,008	%		
	Indium (In)	[7440-74-6]	0,0158	%		
	Iron (Fe)	[7439-89-6]	0,0011	%		
	Lead (Pb)	[7439-92-1]	0,044	%		
	Mercury (Hg)	[7439-97-6]	0,0084	%		
	Nickel (Ni)	[7440-02-0]	0,0159	%		
	Lead (Pb)	[7439-92-1]	0,0047	%		
	Selenium (Se)	[7782-49-2]	~4	ppm		
	Silver (Ag)	[7440-22-4]	4,11	%		
	Zinc (Zn)	[7440-66-6]	~15	ppm		