

Custom Solution - 4 components in Water

Art. ID LM24-CUS-72438-4X100ML
Unit 4 x 100 mL bottle
Deliverydetails No Dangerous Good /not restricted

Description

Conservation 24 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Sodium (Na+)		1000	mg/L		
	Potassium (K+)		1000	mg/L		
	Magnesium (Mg2+)		2500	mg/L		
	Calcium (Ca2+)		4500	mg/L		