

Custom Solution - 4 components 500 mg/L each of in HNO3 5%

Art. ID LM24-CUS-72791
Unit 1 x 100 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation 24 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Niobium (Nb)	[7440-03-1]	500	mg/L		
	Nickel (Ni)	[7440-02-0]	500	mg/L		
	Molybdenum (Mo)	[7439-98-7]	500	mg/L		
	Tin (Sn)	[7440-31-5]	500	mg/L		