

Custom Solution - 5 components 100 µg/mL each of in Acetone

Art. ID LM24-CUS-72794
Unit 1 x 1 mL ampoule
Deliverydetails De Minimis UN1090 II E2

Description

Conservation: 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Imidacloprid D4 (imidazolidin-4,4,5,5 D4)	[1015855-75-0]	100	µg/mL		
	Atrazine D5 (ethylamino D5)	[163165-75-1]	100	µg/mL		
	Terbutylazine D5 (ethyl D5)	[222986-60-9]	100	µg/mL		
	Chlorpyrifos D10 (diethyl D10)	[285138-81-0]	100	µg/mL		
	Isoproturon D6 (dimethyl D6)	[1007461-76-8]	100	µg/mL		