

Custom Solution - 2 components 10 mg/L each of in HNO3 2%

Art. ID LM24-CUS-72797
Unit 1 x 100 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Conservation 12 months (in originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Copper (Cu)	[7440-50-8]	10	mg/L		
	Iron (Fe)	[7439-89-6]	10	mg/L		