

Custom solution - 4 components in Nitric Acid 0.2%

Art. ID LM24-CUS-72807
Unit 100 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG)

Description

Conservation 12 months (originally closed bottle)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Ammonium (NH ₄ ⁺)	[14798-03-9]	100	mg/L		
	Lithium (Li ⁺)		10	mg/L		
	Potassium (K ⁺)		50	mg/L		
	Sodium (Na ⁺)		100	mg/L		