

Custom Solution - Bisphenol A D16 [96210-87-6] 10 mg/L in Methanol

Art. ID LM24-CUS-72827
Unit 3 x 1 mL ampoule
Deliverydetails De Minimis UN1230 (6.1) II E2

Description

Conservation 12 months (originally closed bottle).

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|--------------|---------------------|------|--------|--------|
| | Bisphenol A D16 | [96210-87-6] | 10 | mg/L | | |