

Custom Solution - Bisphenol A D16 [96210-87-6] 1000 µg/L in Methanol + Extended Shelf-Life (Option D): + 12 months

Art. ID LM24-CUS-72829
Unit 10 x 1 mL ampoule
Deliverydetails De Minimis UN1230 (6.1) II E2

Description

Conservation: 12 months + Extended Shelf-Life (Option D) + 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Bisphenol A D16	[96210-87-6]	1000	µg/L		