

Custom Solution - 9 components in HNO3 2%

Art. ID LM24-CUS-72835
Unit 1 x 100 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Conservation 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Chromium (Cr)	[7440-47-3]	50	mg/L		
	Copper (Cu)	[7440-50-8]	50	mg/L		
	Nickel (Ni)	[7440-02-0]	50	mg/L		
	Lead (Pb)	[7439-92-1]	50	mg/L		
	Vanadium (V)	[7440-62-2]	50	mg/L		
	Zinc (Zn)	[7440-66-6]	50	mg/L		
	Arsenic (As)	[7440-38-2]	10	mg/L		
	Cobalt (Co)	[7440-48-4]	10	mg/L		
	Cadmium (Cd)	[7440-43-9]	2	mg/L		