

Custom Solution - 7 components in HNO3 1%, HCl 3% + Extended Shelf-Life (Option D): + 12 months

Art. ID LM24-CUS-72448-250ML
Unit 1 x 250 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Conservation 12 months + Extended Shelf-Life (Option D) + 12 months (originally closed bottle)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Iron (Fe)	[7439-89-6]	500	mg/L		
	Calcium (Ca)	[7440-70-2]	200	mg/L		
	Sulfur (S)	[7704-34-9]	100	mg/L		
	Zinc (Zn)	[7440-66-6]	20	mg/L		
	Manganese (Mn)	[7439-96-5]	20	mg/L		
	Arsenic (As)	[7440-38-2]	1	mg/L		
	Phosphorus (P)	[7723-14-0]	1000	mg/L		