

## **Custom Solution - 25 components 50 mg/L each of in HCl 10%**

Art. ID LM24-CUS-72561-1000ML  
Unit 1 x 1000 mL bottle  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1789 II E2

### Description

Conservation 12 months (in originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	50	mg/L		
	Arsenic (As)	[7440-38-2]	50	mg/L		
	Bismuth (Bi)	[7440-69-9]	50	mg/L		
	Calcium (Ca)	[7440-70-2]	50	mg/L		
	Cadmium (Cd)	[7440-43-9]	50	mg/L		
	Cobalt (Co)	[7440-48-4]	50	mg/L		
	Chromium (Cr)	[7440-47-3]	50	mg/L		
	Copper (Cu)	[7440-50-8]	50	mg/L		
	Iron (Fe)	[7439-89-6]	50	mg/L		
	Potassium (K)	[7440-09-7]	50	mg/L		
	Lithium (Li)	[7439-93-2]	50	mg/L		
	Lanthanum (La)	[7439-91-0]	50	mg/L		
	Magnesium (Mg)	[7439-95-4]	50	mg/L		
	Manganese (Mn)	[7439-96-5]	50	mg/L		
	Molybdenum (Mo)	[7439-98-7]	50	mg/L		
	Sodium (Na)	[7440-23-5]	50	mg/L		
	Nickel (Ni)	[7440-02-0]	50	mg/L		
	Phosphorus (P)	[7723-14-0]	50	mg/L		
	Lead (Pb)	[7439-92-1]	50	mg/L		
	Palladium (Pd)	[7440-05-3]	50	mg/L		
	Rhodium (Rh)	[7440-16-6]	50	mg/L		
	Antimony (Sb)	[7440-36-0]	50	mg/L		
	Strontium (Sr)	[7440-24-6]	50	mg/L		
	Vanadium (V)	[7440-62-2]	50	mg/L		
	Zinc (Zn)	[7440-66-6]	50	mg/L		