

### **Custom Solution - 3 components in Water**

Art. ID LM24-CUS-72865  
Unit 1 x 250 mL bottle  
Deliverydetails No Dangerous Good /not restricted

#### **Description**

Conservation 24 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Phosphorus (P)	[7723-14-0]	1000	mg/L		
	Sulfur (S)	[7704-34-9]	1000	mg/L		
	Silicon (Si)	[7440-21-3]	100	mg/L		