

**Custom Solution - 2 Components 5000 mg/L each of in n-Heptane + Extended Shelf-Life (Option D): + 12 months**

Art. ID LM24-CUS-70443-5X1ML  
Unit 5 x 1 mL ampoule  
Deliverydetails De Minimis UN1206 II E2

Description

Conservation 12 months (originally closed bottle) + Extended Shelf-Life (Option D) + 12 months. (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Mineral oil type A	Mineral Oil [CAS:68334-30-5]	[68334-30-5]	5000	mg/L		
Mineral oil type B	Mineral Oil [CAS:74869-22-0]	[74869-22-0]	5000	mg/L		