

3-Epi-25-Hydroxyvitamin D3 Calibration Solution

Art. ID NIST-2968
 Unit 5 ampoules x 1,2 mL
 Deliverydetails Dry ice shipment UN1170 III E1

Description

This Standard Reference Material (SRM) and its associated certified values are intended primarily for use in calibration of instruments and techniques used for the determination of 3-epi-25-hydroxyvitamin D3. A unit of SRM 2968 consists of five, 2-milliliter ampoules of a single ethanolic 3-Epi-25-Hydroxyvitamin D3 Calibration Solution [3-Epi-25(OH)D3 Calibration Solution], each containing approximately 1.2 mL of solution. /// Sample value(s) - please ask for current certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Mass Fraction	3-epi-25-hydroxyvitamin D3	[73809-05-9]	293,4 ± 13,5	ng/g		
Concentration	3-epi-25-hydroxyvitamin D3	[73809-05-9]	577,0 ± 28,5	nmol/L		