

### **3-Epi-25-Hydroxyvitamin D3 Calibration Solution**

Art. ID	NIST-2968
Unit	5 ampoules x 1,2 mL
Deliverydetails	Dry ice shipment UN1170 III E1

#### Description

This Standard Reference Material (SRM) and its associated certified values are intended primarily for use in calibration of instruments and techniques used for the determination of 3-epi-25-hydroxyvitamin D3. A unit of SRM 2968 consists of five, 2-milliliter ampoules of a single ethanolic 3-Epi-25-Hydroxyvitamin D3 Calibration Solution [3-Epi-25(OH)D3 Calibration Solution], each containing approximately 1.2 mL of solution. /// Sample value(s) - please ask for current certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Mass Fraction	3-epi-25-hydroxyvitamin D3	[73809-05-9]	293,4 ± 13,5	ng/g		
Concentration	3-epi-25-hydroxyvitamin D3	[73809-05-9]	577,0 ± 28,5	nmol/L		