

Anhydrous Milk Fat (triglycerides)

Art. ID ERM-BD519
Unit 5 mL
Deliverydetails No Dangerous Good /not restricted

Description

ERM-BD519 is an anhydrous milk fat. This CRM is available in amber glass ampoules containing 5 mL of anhydrous milk fat, which were sealed under an atmosphere of argon.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
triglyceride	C24		0,078 ± 0,009	g/100 g		
triglyceride	C26		0,338 ± 0,026	g/100 g		
triglyceride	C28		0,76 ± 0,03	g/100 g		
triglyceride	C30		1,45 ± 0,04	g/100 g		
triglyceride	C32		3,03 ± 0,07	g/100 g		
triglyceride	C34		6,6 ± 0,14	g/100 g		
triglyceride	C36		11,25 ± 0,29	g/100 g		
triglyceride	C38		12,46 ± 0,22	g/100 g		
triglyceride	C40		9,91 ± 0,16	g/100 g		
triglyceride	C42		7,4 ± 0,09	g/100 g		
triglyceride	C44		7,1 ± 0,09	g/100 g		
triglyceride	C46		7,74 ± 0,11	g/100 g		
triglyceride	C48		9,07 ± 0,14	g/100 g		
triglyceride	C50		10,34 ± 0,19	g/100 g		
triglyceride	C52		8,33 ± 0,19	g/100 g		
triglyceride	C54		3,83 ± 0,09	g/100 g		
	Cholesterol	[57-88-5]	0,303 ± 0,022	g/100 g		