

## **Vis Wavecal Solution - 8 components (Standard Equivalent to Perkin Elmer)**

Art. ID LM24-OP-4013398  
Unit 500 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

### Description

Conservation 36 months (in originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Barium (Ba)	[7440-39-3]	1	mg/L		
	Calcium (Ca)	[7440-70-2]	1	mg/L		
	Potassium (K)	[7440-09-7]	50	mg/L		
	Lanthanum (La)	[7439-91-0]	10	mg/L		
	Lithium (Li)	[7439-93-2]	10	mg/L		
	Manganese (Mn)	[7439-96-5]	10	mg/L		
	Sodium (Na)	[7440-23-5]	10	mg/L		
	Strontium (Sr)	[7440-24-6]	10	mg/L		