

Silicon Dioxide SiO₂ - 1000 mg/L in H₂O/tr. NaOH for ICP

Art. ID LM24-OP-4013400
Unit 100 mL
Deliverydetails No Dangerous Good /not restricted

Description

Validity: 36 months /// Produced under ISO 17034 and ISO 17025

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------------------------|-------------|---------------------|------|--------|--------|
| | Silicon dioxide (SiO ₂) | [7631-86-9] | 1000 | mg/L | | |