

Custom Solution - 12 components 20 mg/l each of in HNO3 5%

Art. ID LM24-CUS-72929

Unit 1 x 100 mL bottle

Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Bismuth (Bi)	[7440-69-9]	20	mg/L		
	Tellurium (Te)	[13494-80-9]	20	mg/L		
	Niobium (Nb)	[7440-03-1]	20	mg/L		
	Gallium (Ga)	[7440-55-3]	20	mg/L		
	Tantalum (Ta)	[7440-25-7]	20	mg/L		
	Germanium (Ge)	[7440-56-4]	20	mg/L		
	Tin (Sn)	[7440-31-5]	20	mg/L		
	Antimony (Sb)	[7440-36-0]	20	mg/L		
	Tungsten (W)	[7440-33-7]	20	mg/L		
	Zirconium (Zr)	[7440-67-7]	20	mg/L		
	Boron (B)	[7440-42-8]	20	mg/L		
	Titanium (Ti)	[7440-32-6]	20	mg/L		