

### **Custom Solution - 1,2-Dichloroethane D4 [17060-07-0] 5000 mg/L in Methanol**

Art. ID LM24-CUS-72933  
Unit 3 x 1 mL ampoule  
Deliverydetails De Minimis UN1230 (6.1) II E2

#### Description

CRM, ISO 17034 and 17025. Conservation 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	1,2-Dichloroethane D4	[17060-07-0]	5000	mg/L		