

Custom Solution - Phenanthrene D10 [1517-22-2] 1000 mg/L in Dichloromethane

Art. ID LM24-CUS-72934
Unit 3 x 1 mL ampoule
Deliverydetails De Minimis UN1593 III E1

Description

CRM, ISO 17034 and 17025. Conservation 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Phenanthrene D10	[1517-22-2]	1000	mg/L		