

### **Custom Solution - Acenaphthene D10 [15067-26-2] 1000 mg/L in Dichloromethane**

Art. ID LM24-CUS-72935  
Unit 3 x 1 mL ampoule  
Deliverydetails De Minimis UN1593 III E1

#### Description

CRM, ISO 17034 and 17025. Conservation 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Acenaphthene D10	[15067-26-2]	1000	mg/L		