

## **Custom Solution - 31 components 100 mg/L each of in HNO3 5%**

Art. ID LM24-CUS-72936

Unit 1 x 100 mL bottle

Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

### Description

Conservation 24 months (originally closed bottle).

| Text/Information | Analyte/Parameter | CAS number  | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
|                  | Silver (Ag)       | [7440-22-4] | 100                 | mg/L |        |        |
|                  | Aluminium (Al)    | [7429-90-5] | 100                 | mg/L |        |        |
|                  | Arsenic (As)      | [7440-38-2] | 100                 | mg/L |        |        |
|                  | Boron (B)         | [7440-42-8] | 100                 | mg/L |        |        |
|                  | Barium (Ba)       | [7440-39-3] | 100                 | mg/L |        |        |
|                  | Beryllium (Be)    | [7440-41-7] | 100                 | mg/L |        |        |
|                  | Calcium (Ca)      | [7440-70-2] | 100                 | mg/L |        |        |
|                  | Cadmium (Cd)      | [7440-43-9] | 100                 | mg/L |        |        |
|                  | Cobalt (Co)       | [7440-48-4] | 100                 | mg/L |        |        |
|                  | Chromium (Cr)     | [7440-47-3] | 100                 | mg/L |        |        |
|                  | Copper (Cu)       | [7440-50-8] | 100                 | mg/L |        |        |
|                  | Iron (Fe)         | [7439-89-6] | 100                 | mg/L |        |        |
|                  | Lithium (Li)      | [7439-93-2] | 100                 | mg/L |        |        |
|                  | Potassium (K)     | [7440-09-7] | 100                 | mg/L |        |        |
|                  | Magnesium (Mg)    | [7439-95-4] | 100                 | mg/L |        |        |
|                  | Manganese (Mn)    | [7439-96-5] | 100                 | mg/L |        |        |
|                  | Molybdenum (Mo)   | [7439-98-7] | 100                 | mg/L |        |        |
|                  | Sodium (Na)       | [7440-23-5] | 100                 | mg/L |        |        |
|                  | Nickel (Ni)       | [7440-02-0] | 100                 | mg/L |        |        |
|                  | Phosphorus (P)    | [7723-14-0] | 100                 | mg/L |        |        |
|                  | Lead (Pb)         | [7439-92-1] | 100                 | mg/L |        |        |
|                  | Antimony (Sb)     | [7440-36-0] | 100                 | mg/L |        |        |
|                  | Selenium (Se)     | [7782-49-2] | 100                 | mg/L |        |        |
|                  | Tin (Sn)          | [7440-31-5] | 100                 | mg/L |        |        |
|                  | Silicon (Si)      | [7440-21-3] | 100                 | mg/L |        |        |
|                  | Strontium (Sr)    | [7440-24-6] | 100                 | mg/L |        |        |
|                  | Titanium (Ti)     | [7440-32-6] | 100                 | mg/L |        |        |
|                  | Thallium (Tl)     | [7440-28-0] | 100                 | mg/L |        |        |
|                  | Vanadium (V)      | [7440-62-2] | 100                 | mg/L |        |        |

|             |             |     |      |
|-------------|-------------|-----|------|
| Zinc (Zn)   | [7440-66-6] | 100 | mg/L |
| Uranium (U) | [7440-61-1] | 100 | mg/L |