

Custom Solution - 8 components in HNO3 2%

Art. ID LM24-CUS-72938
Unit 1 x 500 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Conservation 8 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Beryllium (Be)	[7440-41-7]	50	µg/L		
	Cadmium (Cd)	[7440-43-9]	50	µg/L		
	Selenium (Se)	[7782-49-2]	100	µg/L		
	Antimony (Sb)	[7440-36-0]	150	µg/L		
	Cobalt (Co)	[7440-48-4]	500	µg/L		
	Arsenic (As)	[7440-38-2]	500	µg/L		
	Lead (Pb)	[7439-92-1]	5000	µg/L		
	Thallium (Tl)	[7440-28-0]	25	µg/L		