

Custom Solution - 7 components 2 mg/L each of in HNO3 2%

Art. ID LM24-CUS-72942
Unit 1 x 100 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Conservation 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Germanium (Ge)	[7440-56-4]	2	mg/L		
	Molybdenum (Mo)	[7439-98-7]	2	mg/L		
	Antimony (Sb)	[7440-36-0]	2	mg/L		
	Silicon (Si)	[7440-21-3]	2	mg/L		
	Tin (Sn)	[7440-31-5]	2	mg/L		
	Titanium (Ti)	[7440-32-6]	2	mg/L		
	Sulfur (S)	[7704-34-9]	2	mg/L		