

Silicon Molybdenum Cast Iron, RM, Block 48 mm x 42 mm x 12 mm

Art. ID BAS-SIMO 2/4
Unit disc
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	0,013 ± 0,002	%		
	Arsenic (As)	[7440-38-2]	0,002 ± 0,001	%		
	Carbon (C)	[7440-44-0]	2,26 ± 0,02	%		
	Cerium (Ce)	[7440-45-1]	0,006 ± 0,001	%		
	Chromium (Cr)	[7440-47-3]	0,896 ± 0,018	%		
	Cobalt (Co)	[7440-48-4]	0,005 ± 0,001	%		
	Copper (Cu)	[7440-50-8]	0,009 ± 0,001	%		
	Manganese (Mn)	[7439-96-5]	0,44 ± 0,005	%		
	Manganese (Mn)	[7439-96-5]	0,002 ± 0,001	%		
	Molybdenum (Mo)	[7439-98-7]	0,443 ± 0,005	%		
	Nickel (Ni)	[7440-02-0]	0,015 ± 0,001	%		
	Phosphorus (P)	[7723-14-0]	0,026 ± 0,001	%		
	Silicon (Si)	[7440-21-3]	4,49 ± 0,01	%		
	Sulfur (S)	[7704-34-9]	0,009 ± 0,001	%		
	Tin (Sn)	[7440-31-5]	0,056 ± 0,002	%		
	Titanium (Ti)	[7440-32-6]	0,006 ± 0,001	%		
	Vanadium (V)	[7440-62-2]	0,008 ± 0,001	%		