

Custom Solution - ICP-OES Wavelength calibration solution - 15 components in HNO3 5%

Art. ID LM24-CUS-70227-2X500ML
Unit 2 x 500 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation 36 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	5	µg/ml		
	Arsenic (As)	[7440-38-2]	5	µg/ml		
	Barium (Ba)	[7440-39-3]	5	µg/ml		
	Cadmium (Cd)	[7440-43-9]	5	µg/ml		
	Cobalt (Co)	[7440-48-4]	5	µg/ml		
	Chromium (Cr)	[7440-47-3]	5	µg/ml		
	Copper (Cu)	[7440-50-8]	5	µg/ml		
	Manganese (Mn)	[7439-96-5]	5	µg/ml		
	Molybdenum (Mo)	[7439-98-7]	5	µg/ml		
	Nickel (Ni)	[7440-02-0]	5	µg/ml		
	Lead (Pb)	[7439-92-1]	5	µg/ml		
	Selenium (Se)	[7782-49-2]	5	µg/ml		
	Strontium (Sr)	[7440-24-6]	5	µg/ml		
	Zinc (Zn)	[7440-66-6]	5	µg/ml		
	Potassium (K)	[7440-09-7]	50	µg/ml		