

AAEM Resin (Acetoacetoxyethyl Ketoester Resin)

Art. ID TRC-A104495-10G
Unit 10 g

Description

Category: Resin, Polymers /// Application Notes: PL-AAEM Resin, is a styrene/methacrylate copolymer matrix. It is a scavenger resin designed to selectively remove the primary amines in the presence of secondary amines. /// References: Yu, Z. R., et al.: Tetrahedron Lett, 41, 8963 (2000)| Yu, Z. R., et al.: J. Chem. Soc, Perkin 1, 1947 (2001)|

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|--|--------------|---------------------|------|--------|--------|
| | AAEM Resin (Acetoacetox yethyl Ketoester Resin) | [34500-18-0] | | | | |