

ABT 925-d6 Fumerate

Art. ID TRC-A112602-50MG
Unit 50 mg

Description

Category: Isotope Labelled Compounds, Inhibitors, Pharmaceuticals, Intermediates and Fine Chemicals, Sulfur and Selenium Compounds /// Application Notes: ABT 925-d6 Fumerate is the isotope labelled analog of ABT 925 Fumerate. ABT 925 Fumerate is a dopamine D3 receptor antagonist that is used in the treatment of acute schizophrenia. /// References: Redden, L., et al.: J. Clin. Psychopharmacol, 31, 221 (2011)|

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	ABT 925-d6 Fumerate					