

Adipoyl-L-carnitine-d3

Art. ID TRC-A291707-25MG
Unit 25 mg

Description

Category: Chiral Reagents, Glycerols, Fatty Acid Derivatives and Lipids, Isotope Labelled Compounds, Metabolites and Impurities /// Application Notes: Adipoyl-L-carnitine-d3 is the labeled analogue of Adipoyl-L-carnitine (A291705), a serum metabolite that varies in concentration in response to normal daily activities including food intake. /// References: Thompson, D.K., et. al.: Metabolomics, 8, 556 (2012)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Adipoyl-L-carnitine-d3					